

Antibodies Product Datasheet

Product Name:	CD294 Antibody
Catalog Number:	C30537
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PTGDR2. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	PTGDR2
Protein Name:	Prostaglandin D2 receptor 2
Other Name:	PTGDR2; CRTH2; DL1R; GPR44; Prostaglandin D2 receptor 2; Chemoattractant receptor-homologous molecule expressed on Th2 cells; G-protein coupled receptor 44; CD294
Human Swiss Prot No.:	Q9Y5Y4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD294 Polyclonal Antibody detects endogenous levels of CD294 protein.
Storage and Stability:	-20°C/1 year