

Antibodies Product Datasheet

Product Name:	Catenin- γ Monoclonal Antibody
Catalog Number:	N0094
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Catenin- γ expressed in E. Coli.
Concentration:	
Gene Name:	JUP
Protein Name:	Junction plakoglobin
Other Name:	JUP; CTNNG; DP3; Junction plakoglobin; Catenin gamma; Desmoplakin III; Desmoplakin-3
Human Swiss Prot No.:	P14923
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Catenin- γ Monoclonal Antibody detects endogenous levels of Catenin- γ protein.
Storage and Stability:	-20°C/1 year