

Antibodies Product Datasheet

Product Name:	CLP36 Antibody
Catalog Number:	C10012
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PDLIM1. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	PDLIM1
Protein Name:	PDZ and LIM domain protein 1
Other Name:	PDLIM1; CLIM1; CLP36; PDZ and LIM domain protein 1; C-terminal LIM domain protein 1; Elfin; LIM domain protein CLP-36
Human Swiss Prot No.:	O00151
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLP36 Polyclonal Antibody detects endogenous levels of CLP36 protein.
Storage and Stability:	-20°C/1 year