

Antibodies Product Datasheet

Product Name:	Neuroglycan C Antibody
Catalog Number:	C15097
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSPG5. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	CSPG5
Protein Name:	Chondroitin sulfate proteoglycan 5
Other Name:	CSPG5; CALEB; NGC; Chondroitin sulfate proteoglycan 5; Acidic leucine-rich EGF-like domain-containing brain protein; Neuroglycan C
Human Swiss Prot No.:	O95196
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Neuroglycan C Polyclonal Antibody detects endogenous levels of Neuroglycan C protein.
Storage and Stability:	-20°C/1 year