

Antibodies Product Datasheet

Product Name:	PI4K2B Antibody
Catalog Number:	C50899
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PI4K2B
Concentration:	1 mg/ml
Gene Name:	PI4K2B
Protein Name:	Phosphatidylinositol 4-kinase type 2-beta (EC 2.7.1.67) (Phosphatidylinositol 4-kinase type II-beta) (PI4KII-BETA)
Other Name:	
Human Swiss Prot No.:	Q8TCG2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PI4K2B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year