

Antibodies Product Datasheet

Product Name:	CDK5 (Phospho-Ser159) Antibody
Catalog Number:	A8496
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CDK5 (Phospho-Ser159)
Concentration:	1 mg/ml
Gene Name:	CDK5 CDKN5
Protein Name:	CDK5 (Phospho-Ser159)
Other Name:	Cyclin-dependent kinase 5 (EC 2.7.11.22) (Cell division protein kinase 5) (Serine/threonine-protein kinase PSSALRE) (Tau protein kinase II catalytic subunit) (TPKII catalytic subunit)
Human Swiss Prot No.:	Q00535
Formulation:	Liquid in PBS containing 50% glycerol, and 0.123% sodium azide.
Specificity:	This antibody detects endogenous levels of CDK5 (Phospho-Ser159) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year