
Antibodies Product Datasheet

Product Name:	C32 Antibody
Catalog Number:	C30733
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Chemokine-like factor at AA range: 103-152
Concentration:	1 mg/ml
Gene Name:	CKLF
Protein Name:	Chemokine-like factor
Other Name:	CKLF; CKLF1; HSPC224; Chemokine-like factor; C32
Human Swiss Prot No.:	Q9UBR5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C32 Polyclonal Antibody detects endogenous levels of C32 protein.
Storage and Stability:	-20°C/1 year