

## Antibodies Product Datasheet

<b>Product Name:</b>	APE1 (Ab-6) Antibody
<b>Catalog Number:</b>	D12044N
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	APE1 (Ab-6) Antibody
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific peptide.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human APE1 around the non-acetylation site of Lys6.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	APEX1
<b>Other Name:</b>	APEX1; APE; APE1; APEX; APX; HAP1; REF1; DNA-(apurinic or apyrimidinic site) lyase; APEX nuclease; APEN; Apurinic-apyrimidinic endonuclease 1; AP endonuclease 1; APE-1; REF-1; Redox factor-1
<b>Human Swiss Prot No.:</b>	P27695
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	APE1 (Ab-6) Antibody detects endogenous levels of APE16 protein.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.