

Antibodies Product Datasheet

Product Name:	CDKN3 Antibody
Catalog Number:	C11137sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CDKN3 Polyclonal
Concentration:	1 mg/ml
Gene Name:	CDKN3 CDI1 CIP2 KAP
Protein Name:	Cyclin-dependent kinase inhibitor 3 (EC 3.1.3.16) (EC 3.1.3.48) (CDK2-associated dual-specificity phosphatase) (Cyclin-dependent kinase interactor 1) (Cyclin-dependent kinase-interacting protein 2) (Kinase-associated phosphatase)
Other Name:	Cyclin-dependent kinase inhibitor 3 (EC 3.1.3.16) (EC 3.1.3.48) (CDK2-associated dual-specificity phosphatase) (Cyclin-dependent kinase interactor 1) (Cyclin-dependent kinase-interacting protein 2) (Kinase-associated phosphatase)
Human Swiss Prot No.:	Q16667
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CDKN3.
Storage and Stability:	-20°C/1 year