

## Antibodies Product Datasheet

<b>Product Name:</b>	RBP4 mouse mAb
<b>Catalog Number:</b>	N1042
<b>Reactivity:</b>	Human
<b>Applications:</b>	ELISA
<b>Dilutions:</b>	ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human RBP4 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	rbp4
<b>Protein Name:</b>	
<b>Other Name:</b>	OTTHUMP00000020114;OTTHUMP00000020115;OTTHUMP00000020116;Plasma retinol binding protein 4;Plasma retinol-binding protein;Plasma retinol-binding protein(1-176);PRBP;PRO2222;RBP;RBP4;RDCCAS;RET4_HUMAN;Retinol binding protein 4;Retinol binding protein 4 inte
<b>Human Swiss Prot No.:</b>	P02753
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	
<b>Storage and Stability:</b>	-20°C/1 year