

Antibodies Product Datasheet

Product Name:	Tryptase-3 Antibody
Catalog Number:	C30867
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Tryptase-3 at AA range: 51-100
Concentration:	1 mg/ml
Gene Name:	TPSD1
Protein Name:	Tryptase-3
Other Name:	Tryptase delta (EC 3.4.21.59) (Delta-tryptase) (HmMCP-3-like tryptase III) (Mast cell mMCP-7-like) (Tryptase-3)
Human Swiss Prot No.:	Q9BZJ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Tryptase-3 Polyclonal Antibody detects endogenous levels of Tryptase-3
Storage and Stability:	-20°C/1 year