

Antibodies Product Datasheet

Product Name:	α -tubulin Monoclonal Antibody(8F11)
Catalog Number:	N1299
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	WB: 1:500-10000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Tubulin alpha-1A chain
Concentration:	
Gene Name:	TUBA1A
Protein Name:	Tubulin alpha-1A chain
Other Name:	
Human Swiss Prot No.:	Q71U36/P68363
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous α -tubulin protein.
Storage and Stability:	-20°C/1 year