

Antibodies Product Datasheet

Product Name:	Smad3 Antibody
Catalog Number:	B1003
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Smad3. AA range:145-194
Concentration:	1 mg/ml
Gene Name:	SMAD3
Protein Name:	Mothers against decapentaplegic homolog 3
Other Name:	SMAD3; MADH3; Mothers against decapentaplegic homolog 3; MAD homolog 3; Mad3; Mothers against DPP homolog 3; hMAD-3; JV15-2; SMAD family member 3; SMAD 3; Smad3; hSMAD3
Human Swiss Prot No.:	P84022
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Smad3 Polyclonal Antibody detects endogenous levels of Smad3 protein.
Storage and Stability:	-20°C/1 year