

Antibodies Product Datasheet

Product Name:	Anti-Human S100A1 (ABT226) Mouse Monoclonal Antibody
Catalog Number:	V0450
Clone ID:	ABT226
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	S100A1
Gene Name:	P23297
Other Name:	Protein S100 A13, Protein S100-A13, S100 calcium binding protein A13, S100 calcium-binding protein A13, S100a13, S10AD_HUMAN
Isotype	IgG1, Kappa
Localization:	Tonsil/ Kindey
Positive Tissue Control:	Nuclear, Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human S100A1 protein, and has no cross reaction with S100A4 or S100A6.
Storage and Stability:	-20°C/1 year