

Antibodies Product Datasheet

Product Name:	Histone deacetylase 10 Antibody
Catalog Number:	C0222
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HDAC10. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	HDAC10
Protein Name:	Histone deacetylase 10
Other Name:	HDAC10; Histone deacetylase 10; HD10
Human Swiss Prot No.:	Q969S8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Histone deacetylase 10 Polyclonal Antibody detects endogenous levels of Histone deacetylase 10 protein.
Storage and Stability:	-20°C/1 year