

Antibodies Product Datasheet

Product Name:	CD240d Antibody
Catalog Number:	C30510
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Blood group Rh(D) polypeptide at AA range: 161-210
Concentration:	1 mg/ml
Gene Name:	RHD
Protein Name:	Blood group Rh(D) polypeptide
Other Name:	RHD; Blood group Rh(D) polypeptide; RHXIII; Rh polypeptide 2; RhPII; Rhesus D antigen; CD240D
Human Swiss Prot No.:	Q02161
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD240d Polyclonal Antibody detects endogenous levels of CD240d protein.
Storage and Stability:	-20°C/1 year