

## Antibodies Product Datasheet

<b>Product Name:</b>	AP-1 Monoclonal Antibody
<b>Catalog Number:</b>	N0027
<b>Reactivity:</b>	Human, Mouse, Monkey
<b>Applications:</b>	WB, IHC-p, IF/ICCFCM, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human AP-1 expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	JUN
<b>Protein Name:</b>	Transcription factor AP-1
<b>Other Name:</b>	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
<b>Human Swiss Prot No.:</b>	P05412
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
<b>Specificity:</b>	AP-1 Monoclonal Antibody detects endogenous levels of AP-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year