

Antibodies Product Datasheet

Product Name:	Cytokeratin 7 Monoclonal Antibody
Catalog Number:	N0188
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cytokeratin 7 expressed in E. Coli.
Concentration:	
Gene Name:	KRT7
Protein Name:	Keratin type II cytoskeletal 7
Other Name:	KRT7; SCL; Keratin; type II cytoskeletal 7; Cytokeratin-7; CK-7; Keratin-7; K7; Sarcolectin; Type-II keratin Kb7
Human Swiss Prot No.:	P08729
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cytokeratin 7 Monoclonal Antibody detects endogenous levels of Cytokeratin 7 protein.
Storage and Stability:	-20°C/1 year