

Antibodies Product Datasheet

Product Name:	Flg Monoclonal Antibody
Catalog Number:	N0275
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human Flg (aa22-376) fused with hlgGFc tag expressed in HEK293 cells.
Concentration:	
Gene Name:	FGFR1
Protein Name:	Basic fibroblast growth factor receptor 1
Other Name:	FGFR1; BFGFR; CEK; FGFBR; FLG; FLT2; HBGFR; Fibroblast growth factor receptor 1; FGFR-1; Basic fibroblast growth factor receptor 1; BFGFR; bFGF-R-1; Fms-like tyrosine kinase 2; FLT-2; N-sam; Proto-oncogene c-Fgr; CD antigen CD331
Human Swiss Prot No.:	P11362
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Flg Monoclonal Antibody detects endogenous levels of Flg protein.
Storage and Stability:	-20°C/1 year