
Antibodies Product Datasheet

Product Name:	TGFβ1 mouse Monoclonal Antibody(10E5)
Catalog Number:	N1609
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of TGFβ1 at AA range of 310-390
Concentration:	1 mg/ml
Gene Name:	TGFB1 TGFB
Protein Name:	Transforming growth factor beta-1, TGF-β1, TGF b
Other Name:	TGFB1
Human Swiss Prot No.:	P01137
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TGFβ1 protein detects endogenous levels of TGFB1
Storage and Stability:	-20°C/1 year