
Antibodies Product Datasheet

Product Name:	Maf Antibody
Catalog Number:	R12-2985
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Maf Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Maf.
Concentration:	1 mg/ml
Gene Name:	MAF
Other Name:	Transcription factor Maf; Proto-oncogene c-maf; C-Maf; V-maf musculoaponeurotic fibrosarcoma oncogene homolog; MAF
Human Swiss Prot No.:	O75444
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Maf Antibodyantibody detects endogenous levels of Maf.
Storage and Stability:	Stable at -20°C for at least 1 year.