

Antibodies Product Datasheet

Product Name:	CLMP Antibody
Catalog Number:	C49893
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CLMP
Concentration:	1 mg/ml
Gene Name:	CLMP ACAM ASAM UNQ318/PRO363
Protein Name:	CXADR-like membrane protein (Adipocyte adhesion molecule) (Coxsackie- and adenovirus receptor-like membrane protein) (CAR-like membrane protein)
Other Name:	
Human Swiss Prot No.:	Q9H6B4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CLMP at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year