

Antibodies Product Datasheet

Product Name:	Bcl-10 Monoclonal Antibody
Catalog Number:	N0056
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Bcl-10 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	BCL10
Protein Name:	B-cell lymphoma/leukemia 10
Other Name:	BCL10; CIPER; CLAP; B-cell lymphoma/leukemia 10; B-cell CLL/lymphoma 10; Bcl-10; CARD-containing molecule enhancing NF-kappa-B; CARD-like apoptotic protein; hCLAP; CED-3/ICH-1 prodomain homologous E10-like regulator; CIPER; Cellular homolog
Human Swiss Prot No.:	O95999
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	Bcl-10 Monoclonal Antibody detects endogenous levels of Bcl-10 protein.
Storage and Stability:	-20°C/1 year