

Antibodies Product Datasheet

Product Name:	Fibulin-5 Antibody
Catalog Number:	C15754
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FBLN5. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	FBLN5
Protein Name:	Fibulin-5
Other Name:	FBLN5; DANCE; Fibulin-5; FIBL-5; Developmental arteries and neural crest EGF-like protein; Dance; Urine p50 protein; UP50
Human Swiss Prot No.:	Q9UBX5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fibulin-5 Polyclonal Antibody detects endogenous levels of Fibulin-5 protein.
Storage and Stability:	-20°C/1 year