

Antibodies Product Datasheet

Product Name:	ZNF608 Antibody
Catalog Number:	C19605
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF608. AA range:1231-1280
Concentration:	1 mg/ml
Gene Name:	ZNF608
Protein Name:	Zinc finger protein 608
Other Name:	ZNF608; KIAA1281; Zinc finger protein 608; Renal carcinoma antigen NY-REN-36
Human Swiss Prot No.:	Q9ULD9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ZNF608 Polyclonal Antibody detects endogenous levels of ZNF608 protein.
Storage and Stability:	-20°C/1 year