
Antibodies Product Datasheet

Product Name:	FGFR2 (Phospho-Tyr769) Antibody
Catalog Number:	A8102
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGFR2 (Phospho-Tyr769)
Concentration:	1 mg/ml
Gene Name:	FGFR2 BEK KGFR KSAM
Protein Name:	FGFR2 (Phospho-Tyr769)
Other Name:	Fibroblast growth factor receptor 2 (FGFR-2) (EC 2.7.10.1) (K-sam) (KGFR) (Keratinocyte growth factor receptor) (CD antigen CD332)
Human Swiss Prot No.:	P21802
Formulation:	Liquid in PBS containing 50% glycerol, and 0.128% sodium azide.
Specificity:	This antibody detects endogenous levels of FGFR2 (Phospho-Tyr769) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year