

Antibodies Product Datasheet

Product Name:	DMPK Antibody
Catalog Number:	R12-2704
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	DMPK Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DMPK.
Concentration:	1 mg/ml
Gene Name:	DMPK
Other Name:	Myotonin-protein kinase; Myotonic dystrophy protein kinase; Myotonic dystrophy protein kinase; MDPK; DM-kinase; DMK; DMPK; MT-PK; DMPK; MDPK
Human Swiss Prot No.:	Q09013
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	DMPK Antibodyantibody detects endogenous levels of DMPK.
Storage and Stability:	Stable at -20°C for at least 1 year.