
Antibodies Product Datasheet

Product Name:	Calreticulin Rabbit Antibody
Catalog Number:	C42106
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein derived from human Calreticulin
Concentration:	1 mg/ml
Gene Name:	CALR CRTC
Protein Name:	Calreticulin
Other Name:	Calreticulin (CRP55) (Calregulin) (Endoplasmic reticulum resident protein 60) (ERp60) (HACBP) (grp60); YM3802
Human Swiss Prot No.:	P27797
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Calreticulin
Storage and Stability:	-20°C/1 year