

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Acetyl-Lys proteins Antibody

Catalog Number: D0020

Reactivity: Human, Mouse, Rat, Monkey **Applications:** WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in

other applications.

Modification: Acetyl Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: Synthesized acetyl-peptide derived from human acetylation Lys proteins.

Concentration: 1 mg/ml

Gene Name:

Protein Name:

Other Name:

Human Swiss Prot No.:

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Acetyl-Lys proteins Polyclonal Antibody detects endogenous levels of

acetylated Lys proteins.

Storage and Stability: -20°C/1 year