

Antibodies Product Datasheet

Product Name:	ATGL Antibody
Catalog Number:	C48048
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ATGL
Concentration:	1 mg/ml
Gene Name:	PNPLA2 ATGL FP17548
Protein Name:	Patatin-like phospholipase domain-containing protein 2 (EC 3.1.1.3) (Adipose triglyceride lipase) (Calcium-independent phospholipase A2) (Desnutrin) (IPLA2-zeta) (Pigment epithelium-derived factor) (TTS2.2) (Transport-secretion protein 2) (TTS2)
Other Name:	
Human Swiss Prot No.:	Q96AD5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATGL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year