

## Antibodies Product Datasheet

<b>Product Name:</b>	PR (Phospho-Ser400) Antibody
<b>Catalog Number:</b>	A0559
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Progesterone Receptor around the phosphorylation site of Ser400. AA range:371-420
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PGR
<b>Protein Name:</b>	Progesterone receptor
<b>Other Name:</b>	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
<b>Human Swiss Prot No.:</b>	P06401
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-PR (S400) Polyclonal Antibody detects endogenous levels of PR protein only when phosphorylated at S400.
<b>Storage and Stability:</b>	-20°C/1 year