

Antibodies Product Datasheet

Product Name:	MOB1 (Phospho-Thr12) Antibody
Catalog Number:	A9588
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human MOB1 (Thr12)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	MOB1 (Thr12)
Other Name:	MOB kinase activator 1A (Mob1 alpha) (Mob1A) (Mob1 homolog 1B) (Mps one binder kinase activator-like 1B)
Human Swiss Prot No.:	Q9H8S9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat MOB1 (phospho-Thr12)
Storage and Stability:	-20°C/1 year