

Antibodies Product Datasheet

Product Name:	EMR2 (Cleaved-Leu517) Antibody
Catalog Number:	L0407
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EMR2 (Cleaved-Leu517)
Concentration:	1 mg/ml
Gene Name:	EMR2
Protein Name:	EMR2 (Cleaved-Leu517)
Other Name:	EGF-like module-containing mucin-like hormone receptor-like 2 (EGF-like module receptor 2) (CD antigen CD312)
Human Swiss Prot No.:	Q9UHX3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human EMR2 (Cleaved-Leu517, protein was cleaved amino acid sequence between 517-518)
Storage and Stability:	-20°C/1 year