

## Antibodies Product Datasheet

<b>Product Name:</b>	ABCB9 Antibody
<b>Catalog Number:</b>	C42428
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ABCB9
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ABCB9 KIAA1520
<b>Protein Name:</b>	ABCB9
<b>Other Name:</b>	ATP-binding cassette sub-family B member 9 (ATP-binding cassette transporter 9) (ABC transporter 9 protein) (hABCB9) (TAP-like protein) (TAPL)
<b>Human Swiss Prot No.:</b>	Q9NP78
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.312% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of ABCB9 at Human/Mouse/Rat
<b>Storage and Stability:</b>	-20°C/1 year