
Antibodies Product Datasheet

Product Name:	CD33 Monoclonal Antibody
Catalog Number:	N0118
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD33 (48-258) expressed in E. Coli.
Concentration:	
Gene Name:	CD33
Protein Name:	Myeloid cell surface antigen CD33
Other Name:	CD33; SIGLEC3; Myeloid cell surface antigen CD33; Sialic acid-binding Ig-like lectin 3; Siglec-3; gp67; CD antigen CD33
Human Swiss Prot No.:	P20138
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD33 Monoclonal Antibody detects endogenous levels of CD33 protein.
Storage and Stability:	-20°C/1 year