

## Antibodies Product Datasheet

<b>Product Name:</b>	MINK1 Antibody
<b>Catalog Number:</b>	C11342
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human MAP4K6. AA range:401-450
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MINK1
<b>Protein Name:</b>	Misshapen-like kinase 1
<b>Other Name:</b>	MINK1; B55; MAP4K6; MINK; YSK2; ZC3; Misshapen-like kinase 1; GCK family kinase MiNK; MAPK/ERK kinase kinase 6; MEK kinase kinase 6; MEKKK 6; Misshapen/NIK-related kinase; Mitogen-activated protein kinase kinase kinase 6
<b>Human Swiss Prot No.:</b>	Q8N4C8
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	MINK1 Polyclonal Antibody detects endogenous levels of MINK1 protein.
<b>Storage and Stability:</b>	-20°C/1 year