
Antibodies Product Datasheet

Product Name:	IL-2 Antibody
Catalog Number:	R12-2205
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-2. AA range:16-65
Concentration:	1 mg/ml
Gene Name:	IL2
Protein Name:	Interleukin-2
Other Name:	IL2; Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin
Human Swiss Prot No.:	P60568
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-2 Polyclonal Antibody detects endogenous levels of IL-2 protein.
Storage and Stability:	-20°C/1 year