

Antibodies Product Datasheet

Product Name:	GGT5 (heavy chain, Cleaved-Gly387) Antibody
Catalog Number:	L0430
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GGT5 (heavy chain, Cleaved-Gly387)
Concentration:	1 mg/ml
Gene Name:	GGT5 GGTLA1
Protein Name:	GGT5 (heavy chain, Cleaved-Gly387)
Other Name:	Gamma-glutamyltransferase 5 (GGT 5) (EC 2.3.2.2) (Gamma-glutamyl transpeptidase-related enzyme) (GGT-rel) (Gamma-glutamyltransferase-like activity 1) (Gamma-glutamyltranspeptidase 5) (Glutathione hydrolase 5) (EC 3.4.19.13) (Leukotriene-C4 hydrolase) (EC P36269
Human Swiss Prot No.:	P36269
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat GGT5 (heavy chain, Cleaved-Gly387, protein was cleaved amino acid sequence between 387-388)
Storage and Stability:	-20°C/1 year