
Antibodies Product Datasheet

Product Name:	CD73 Monoclonal Antibody
Catalog Number:	N0132
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD73 expressed in E. Coli.
Concentration:	
Gene Name:	NT5E
Protein Name:	5'-nucleotidase
Other Name:	NT5E; NT5; NTE; 5'-nucleotidase; 5'-NT; Ecto-5'-nucleotidase; CD antigen CD73
Human Swiss Prot No.:	P21589
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD73 Monoclonal Antibody detects endogenous levels of CD73 protein.
Storage and Stability:	-20°C/1 year