
Antibodies Product Datasheet

Product Name:	P53 mouse mAb(PT0337)
Catalog Number:	N1731
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),WB
Dilutions:	IHC-p 1:100-300, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human P53
Concentration:	1 mg/ml
Gene Name:	TP53 P53
Protein Name:	P53
Other Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human P53
Storage and Stability:	-20°C/1 year