

Antibodies Product Datasheet

Product Name:	FADS1 Antibody
Catalog Number:	C51069
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FADS1
Concentration:	1 mg/ml
Gene Name:	FADS1 FADSD5
Protein Name:	Fatty acid desaturase 1 (EC 1.14.19.-) (Delta(5) fatty acid desaturase) (D5D) (Delta(5) desaturase) (Delta-5 desaturase)
Other Name:	
Human Swiss Prot No.:	O60427
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FADS1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year