

Antibodies Product Datasheet

Product Name:	PRDM3 Antibody
Catalog Number:	C10980
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MECOM. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	MECOM
Protein Name:	MDS1 and EVI1 complex locus protein MDS1
Other Name:	MECOM; MDS1; MDS1 and EVI1 complex locus protein MDS1; Myelodysplasia syndrome 1 protein; Myelodysplasia syndrome-associated protein 1
Human Swiss Prot No.:	Q13465
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PRDM3 Polyclonal Antibody detects endogenous levels of PRDM3 protein.
Storage and Stability:	-20°C/1 year