
Antibodies Product Datasheet

Product Name:	AR Monoclonal Antibody
Catalog Number:	N0042
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of AR (aa689-919) expressed in E. Coli.
Concentration:	
Gene Name:	AR
Protein Name:	Androgen receptor
Other Name:	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Human Swiss Prot No.:	P10275
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	AR Monoclonal Antibody detects endogenous levels of AR protein.
Storage and Stability:	-20°C/1 year