

Antibodies Product Datasheet

Product Name:	Cdk4 Antibody
Catalog Number:	C30140
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CDK4. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CDK4
Protein Name:	Cyclin-dependent kinase 4
Other Name:	CDK4; Cyclin-dependent kinase 4; Cell division protein kinase 4; PSK-J3
Human Swiss Prot No.:	P11802
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdk4 Polyclonal Antibody detects endogenous levels of Cdk4 protein.
Storage and Stability:	-20°C/1 year