
Antibodies Product Datasheet

Product Name:	Sin1 (Phospho-Thr86) Antibody
Catalog Number:	A9623
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Sin1 (Thr86)
Concentration:	1 mg/ml
Gene Name:	MAPKAP1 MIP1 SIN1
Protein Name:	Sin1 (Thr86)
Other Name:	Target of rapamycin complex 2 subunit MAPKAP1 (TORC2 subunit MAPKAP1) (Mitogen-activated protein kinase 2-associated protein 1) (Stress-activated map kinase-interacting protein 1) (SAPK-interacting protein 1) (mSIN1)
Human Swiss Prot No.:	Q9BPZ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Sin1 (phospho-Thr86)
Storage and Stability:	-20°C/1 year