

Antibodies Product Datasheet

Product Name:	HCP1/PCFT Antibody
Catalog Number:	C51067
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HCP1/PCFT
Concentration:	1 mg/ml
Gene Name:	SLC46A1 HCP1 PCFT
Protein Name:	Proton-coupled folate transporter (G21) (Heme carrier protein 1) (PCFT/HCP1) (Solute carrier family 46 member 1)
Other Name:	
Human Swiss Prot No.:	Q96NT5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HCP1/PCFT at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year