
Antibodies Product Datasheet

Product Name:	p16 INK mouse mAb(PTR2302)
Catalog Number:	N1754
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Mouse, Monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p16 INK
Concentration:	
Gene Name:	CDKN2A CDKN2 MTS1
Protein Name:	p16 INK
Other Name:	Cyclin-dependent kinase inhibitor 2A, isoforms 1/2/3 (Cyclin-dependent kinase 4 inhibitor A) (CDK4I) (Multiple tumor suppressor 1) (MTS-1) (p16-INK4a) (p16-INK4) (p16INK4A)
Human Swiss Prot No.:	P42771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.140% sodium azide.
Specificity:	This antibody detects endogenous levels of p16 INK at Human
Storage and Stability:	-20°C/1 year