

Antibodies Product Datasheet

Product Name:	Anti-Human Mast Cell Tryptase (ABT194) Mouse Monoclonal Antibody
Catalog Number:	V0390
Clone ID:	ABT194
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	TPSAB1
Gene Name:	Q15661
Other Name:	alpha II, Lung tryptase, Mast cell alpha II tryptase, Mast cell beta I tryptase, Mast cell protease 7, Mast cell protease II, MCP 7, Pituitary tryptase, Skin tryptase, TPS 1, TPS1, TPS2, TPSAB1, TPSAB1 protein, TPSB1, Tryptase 1, Tryptase alpha 1, tryptase alpha I included, Tryptase alpha II, tryptase alpha II included, tryptase alpha included, tryptase alpha/beta 1, Tryptase beta 1, tryptase beta I included, Tryptase I, tryptase I included, Tryptase III, Tryptase skin
Isotype	IgG1, Kappa
Localization:	Tonsil
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Mast cell tryptase protein.
Storage and Stability:	-20°C/1 year