

Antibodies Product Datasheet

Product Name:	CCDC64B Antibody
Catalog Number:	C50787
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCDC64B
Concentration:	1 mg/ml
Gene Name:	CCDC64B BICDR2
Protein Name:	Bicaudal D-related protein 2 (BICD-related protein 2) (BICDR-2) (Coiled-coil domain-containing protein 64B)
Other Name:	
Human Swiss Prot No.:	A1A5D9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCDC64B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year