

Antibodies Product Datasheet

Product Name:	JAK2/3 (Phospho-Tyr966/939) Antibody
Catalog Number:	A8120
Reactivity:	Human:Y966/939, Mouse:Y966/935, Rat:Y966/935
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized Phospho peptide derived from human JAK2/3.at AA range: T966
Concentration:	1 mg/ml
Gene Name:	JAK2 JAK3
Protein Name:	JAK2/3
Other Name:	Tyrosine-protein kinase JAK2/JAK3 (EC 2.7.10.2) (Janus kinase 2/Janus kinase 3) (JAK-2/JAK-3)
Human Swiss Prot No.:	O60674/P52333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of JAK2/3 (Phospho-Tyr966/939), It doesn't react with total protein.
Storage and Stability:	-20°C/1 year